

**STORE AND
TRANSPORT
CHEMICALS
WITH CONFIDENCE**

**CARBOYS
AND
JERRICANS**

Carboys and Jerricans

Accumax Carboys and Jerricans, meticulously crafted from premium virgin-grade resins to redefine excellence in liquid containment. Designed to surpass your expectations, their superior construction ensures unmatched durability and reliability. These versatile containers cater to a myriad of applications, from storage and mixing to dispensing and transportation of diverse liquids and reagents. Certified leakproof for worry-free transit, they provide a secure shield for your valuable products. Engineered with user convenience in mind, these containers boast ergonomic handles, wide-mouth openings, and a convenient spigot for seamless filling and dispensing.

Available in a range of raw materials and capacities, shapes, and volumes, Accumax Carboys and Jerricans offer a tailored solution for every need. Whether you seek a dependable storage option for sensitive materials or a versatile vessel for mixing and transporting substances, trust in Accumax for unparalleled quality and reliability.



Product Features



Made from US FDA
21 CFR compliant
Raw Material



Individually printed Laser
Engraved graduation with
Lot Number marking for
visual determination.*



Leakproof PP closure
design with Silicone
gasket

Convenient
Integrated
shoulder Handle
for ease of carrying



Anti-Wobbling Mouth to
prevent splashing of
liquids in Jerricans



Rectangular Jerrycan
for efficient use of
storage space



Product Features



Round bottom carboys for homogeneous stirring



Stable base and bottoms with curved inner corners for easy cleaning

Threads on bottles and closures have continuous, straight-shouldered semi-buttress threads.



Unique serial number for easy identification based on customer requirement



Variants

Carboys



10, 15, 20, 25 & 50 litres. (With and without Spigot)

Jerricans



10 litres. and 20 litres. (With and without Anti-Wobbling mouth)

Material



PP/HDPE/LDPE with handle

Ordering Information

Carboy CAT Codes			
CAT Code	Description	Quantity/Pack	Quantity/Case
PP-Autoclavable			
ACB-PP-10	10L Polypropylene Carboy W/o Cap	1	1
ACB-PP-15	15L Polypropylene Carboy W/o Cap	1	1
ACB-PP-20	20L Polypropylene Carboy W/o Cap	1	1
ACB-PP-25	25L Polypropylene Carboy W/o Cap	1	1
ACB-PP-50	50L Polypropylene Carboy W/oCap	1	1
ACB-PP-10-C	10L Polypropylene Carboy W/Cap	1	1
ACB-PP-15-C	15L Polypropylene Carboy W/Cap	1	1
ACB-PP-20-C	20L Polypropylene Carboy W/Cap	1	1
ACB-PP-25-C	25L Polypropylene Carboy W/Cap	1	1
ACB-PP-50-C	50L Polypropylene Carboy W/Cap	1	1
ACB-PP-10-SP	10L Polypropylene Carboy with Spigot W/o Cap	1	1
ACB-PP-15-SP	15L Polypropylene Carboy with Spigot W/o Cap	1	1
ACB-PP-20-SP	20L Polypropylene Carboy with Spigot W/o Cap	1	1
ACB-PP-25-SP	25L Polypropylene Carboy with Spigot W/o Cap	1	1
ACB-PP-50-SP	50L Polypropylene Carboy with Spigot W/o Cap	1	1
ACB-PP-10-SPC	10L Polypropylene Carboy with Spigot W/Cap	1	1
ACB-PP-15-SPC	15L Polypropylene Carboy with Spigot W/Cap	1	1
ACB-PP-20-SPC	20L Polypropylene Carboy with Spigot W/Cap	1	1
ACB-PP-25-SPC	25L Polypropylene Carboy with Spigot W/Cap	1	1
ACB-PP-50-SPC	50L Polypropylene Carboy with Spigot W/Cap	1	1

Ordering Information

Carboy CAT Codes			
CAT Code	Description	Quantity/Pack	Quantity/Case
PP - Sterile			
ACB-PP-10-S	10L Polypropylene Carboy W/o Cap Sterile	1	1
ACB-PP-15-S	15L Polypropylene Carboy W/o Cap Sterile	1	1
ACB-PP-20-S	20L Polypropylene Carboy W/o Cap Sterile	1	1
ACB-PP-25-S	25L Polypropylene Carboy W/o Cap Sterile	1	1
ACB-PP-50-S	50L Polypropylene Carboy W/oCap Sterile	1	1
ACB-PP-10-C-S	10L Polypropylene Carboy W/Cap Sterile	1	1
ACB-PP-15-C-S	15L Polypropylene Carboy W/Cap Sterile	1	1
ACB-PP-20-C-S	20L Polypropylene Carboy W/Cap Sterile	1	1
ACB-PP-25-C-S	25L Polypropylene Carboy W/Cap Sterile	1	1
ACB-PP-50-C-S	50L Polypropylene Carboy W/Cap Sterile	1	1
ACB-PP-10-SPC-S	10L Polypropylene Carboy with Spigot W/o Cap Sterile	1	1
ACB-PP-15-SPC-S	15L Polypropylene Carboy with Spigot W/o Cap Sterile	1	1
ACB-PP-20-SPC-S	20L Polypropylene Carboy with Spigot W/o Cap Sterile	1	1
ACB-PP-25-SPC-S	25L Polypropylene Carboy with Spigot W/o Cap Sterile	1	1
ACB-PP-50-SPC-S	50L Polypropylene Carboy with Spigot W/o Cap Sterile	1	1

Ordering Information

Carboy CAT Codes			
CAT Code	Description	Quantity/Pack	Quantity/Case
HDPE			
ACB-HD-10	10L Polyethylene Carboy W/o Cap	1	1
ACB-HD-15	15L Polyethylene Carboy W/o Cap	1	1
ACB-HD-20	20L Polyethylene Carboy W/o Cap	1	1
ACB-HD-25	25L Polyethylene Carboy W/o Cap	1	1
ACB-HD-50	50L Polyethylene Carboy W/o Cap	1	1
ACB-HD-10-C	10L Polyethylene Carboy W/Cap	1	1
ACB-HD-15-C	15L Polyethylene Carboy W/Cap	1	1
ACB-HD-20-C	20L Polyethylene Carboy W/Cap	1	1
ACB-HD-25-C	25L Polyethylene Carboy W/Cap	1	1
ACB-HD-50-C	50L Polyethylene Carboy W/Cap	1	1
ACB-HD-10-SP	10L Polyethylene Carboy with Spigot W/o Cap	1	1
ACB-HD-15-SP	15L Polyethylene Carboy with Spigot W/o Cap	1	1
ACB-HD-20-SP	20L Polyethylene Carboy with Spigot W/o Cap	1	1
ACB-HD-25-SP	25L Polyethylene Carboy with Spigot W/o Cap	1	1
ACB-HD-50-SP	50L Polyethylene Carboy with Spigot W/o Cap	1	1
ACB-HD-10-SPC	10L Polyethylene Carboy with Spigot W/Cap	1	1
ACB-HD-15-SPC	15L Polyethylene Carboy with Spigot W/Cap	1	1
ACB-HD-20-SPC	20L Polyethylene Carboy with Spigot W/Cap	1	1
ACB-HD-25-SPC	25L Polyethylene Carboy with Spigot W/Cap	1	1
ACB-HD-50-SPC	50L Polyethylene Carboy with Spigot W/Cap	1	1

Ordering Information

Carboy CAT Codes			
CAT Code	Description	Quantity/Pack	Quantity/Case
HDPE - Sterile			
ACB-HD-10-S	10L Polyethylene Carboy W/o Cap, Sterile	1	1
ACB-HD-15-S	15L Polyethylene Carboy W/o Cap, Sterile	1	1
ACB-HD-20-S	20L Polyethylene Carboy W/o Cap, Sterile	1	1
ACB-HD-25-S	25L Polyethylene Carboy W/o Cap, Sterile	1	1
ACB-HD-50-S	50L Polyethylene Carboy W/o Cap, Sterile	1	1
ACB-HD-10-C-S	10L Polyethylene Carboy W/Cap, Sterile	1	1
ACB-HD-15-C-S	15L Polyethylene Carboy W/Cap, Sterile	1	1
ACB-HD-20-C-S	20L Polyethylene Carboy W/Cap, Sterile	1	1
ACB-HD-25-C-S	25L Polyethylene Carboy W/Cap, Sterile	1	1
ACB-HD-50-C-S	50L Polyethylene Carboy W/Cap, Sterile	1	1
ACB-HD-10-SPC-S	10L Polyethylene Carboy with Spigot W/Cap, Sterile	1	1
ACB-HD-15-SPC-S	15L Polyethylene Carboy with Spigot W/Cap, Sterile	1	1
ACB-HD-20-SPC-S	20L Polyethylene Carboy with Spigot W/Cap, Sterile	1	1
ACB-HD-25-SPC-S	25L Polyethylene Carboy with Spigot W/Cap, Sterile	1	1
ACB-HD-50-SPC-S	50L Polyethylene Carboy with Spigot W/Cap, Sterile	1	1

Ordering Information

Carboy CAT Codes			
CAT Code	Description	Quantity/Pack	Quantity/Case
LDPE			
ACB-LD-10	10L Polyethylene Carboy W/o Cap	1	1
ACB-LD-15	15L Polyethylene Carboy W/o Cap	1	1
ACB-LD-20	20L Polyethylene Carboy W/o Cap	1	1
ACB-LD-25	25L Polyethylene Carboy W/o Cap	1	1
ACB-LD-50	50L Polyethylene Carboy W/o Cap	1	1
ACB-LD-10-C	10L Polyethylene Carboy W/Cap	1	1
ACB-LD-15-C	15L Polyethylene Carboy W/Cap	1	1
ACB-LD-20-C	20L Polyethylene Carboy W/Cap	1	1
ACB-LD-25-C	25L Polyethylene Carboy W/Cap	1	1
ACB-LD-50-C	50L Polyethylene Carboy W/Cap	1	1
ACB-LD-10-SP	10L Polyethylene Carboy with Spigot W/o Cap	1	1
ACB-LD-15-SP	15L Polyethylene Carboy with Spigot W/o Cap	1	1
ACB-LD-20-SP	20L Polyethylene Carboy with Spigot W/o Cap	1	1
ACB-LD-25-SP	25L Polyethylene Carboy with Spigot W/o Cap	1	1
ACB-LD-50-SP	50L Polyethylene Carboy with Spigot W/o Cap	1	1
ACB-LD-10-SPC	10L Polyethylene Carboy with Spigot W/Cap	1	1
ACB-LD-15-SPC	15L Polyethylene Carboy with Spigot W/Cap	1	1
ACB-LD-20-SPC	20L Polyethylene Carboy with Spigot W/Cap	1	1
ACB-LD-25-SPC	25L Polyethylene Carboy with Spigot W/Cap	1	1
ACB-LD-50-SPC	50L Polyethylene Carboy with Spigot W/Cap	1	1

Ordering Information

Carboy CAT Codes			
CAT Code	Description	Quantity/Pack	Quantity/Case
LDPE - Sterile			
ACB-LD-10-S	10L Polyethylene Carboy W/o Cap, Sterile	1	1
ACB-LD-15-S	15L Polyethylene Carboy W/o Cap, Sterile	1	1
ACB-LD-20-S	20L Polyethylene Carboy W/o Cap, Sterile	1	1
ACB-LD-25-S	25L Polyethylene Carboy W/o Cap, Sterile	1	1
ACB-LD-50-S	50L Polyethylene Carboy W/o Cap, Sterile	1	1
ACB-LD-10-C-S	10L Polyethylene Carboy W/Cap, Sterile	1	1
ACB-LD-15-C-S	15L Polyethylene Carboy W/Cap, Sterile	1	1
ACB-LD-20-C-S	20L Polyethylene Carboy W/Cap, Sterile	1	1
ACB-LD-25-C-S	25L Polyethylene Carboy W/Cap, Sterile	1	1
ACB-LD-50-C-S	50L Polyethylene Carboy W/Cap, Sterile	1	1
ACB-LD-10-SPC-S	10L Polyethylene Carboy with Spigot W/Cap, Sterile	1	1
ACB-LD-15-SPC-S	15L Polyethylene Carboy with Spigot W/Cap, Sterile	1	1
ACB-LD-20-SPC-S	20L Polyethylene Carboy with Spigot W/Cap, Sterile	1	1
ACB-LD-25-SPC-S	25L Polyethylene Carboy with Spigot W/Cap, Sterile	1	1
ACB-LD-50-SPC-S	50L Polyethylene Carboy with Spigot W/Cap, Sterile	1	1
Caps			
A-83B-C	83B PP Screw Closure with Silicon Gasket	1	1
A-83B-C-B	83B PP Screw Closure with Silicon Gasket, Bulk	5	5
A-83B-C-S	83B PP Screw Closure with Silicon Gasket, Sterile	1	1
A-83B-C-B-S	83B PP Screw Closure with Silicon Gasket, Bulk, Sterile	5	5

Ordering Information

Jerricans CAT Codes			
CAT Code	Description	Quantity/Pack	Quantity/Case
AJB-HD-10	10ltrs Jerrican with Cap	1	1
AJB-HD-10-AW	10ltrs Jerrican with Cap, Anti-Wobbling	1	1
AJB-HD-10-B	10ltrs Jerrican with Cap, Bulk	1	6
AJB-HD-10-AW-B	10ltrs Jerrican with Cap, Anti-Wobbling, Bulk	1	6
AJB-HD-20	20ltrs Jerrican with Cap	1	1
AJB-HD-20-AW	20ltrs Jerrican with Cap, Anti-Wobbling	1	1
AJB-HD-20-B	20ltrs Jerrican with Cap, Bulk	1	4
AJB-HD-20-AW-B	20ltrs Jerrican with Cap, Anti-Wobbling, Bulk	1	4

ACCUMAX

Accumax Lab Devices Pvt. Ltd.
www.accumaxlab.com

ISO 9001:2015
ISO 13485:2016

*We also execute OEM projects. Specifications can be changed without notice for quality improvement.